


# Automatic Analysis of Dual-Channel Droplet Digital PCR Experiments to Detect BRAF-V600 Mutations

Dean Attali

<http://deanattali.com>

Jennifer Bryan Lab @ MSL, UBC




MSc exit seminar  
April 15, 2016



a place of mind  
THE UNIVERSITY OF BRITISH COLUMBIA

# Summer 2014

Jenny, can  
you help  
automate our  
analysis?




Dean, it's  
for you!



# *BRAF* Gene / MAPK Pathway

- B-Raf protein kinase
- Normal conditions:  
**Growth factor** binds  $\Rightarrow$   
**Ras** protein activated  $\Rightarrow$   
**B-Raf** protein activated  $\Rightarrow$   
More phosphorylations  $\Rightarrow$   
Signal for **cell** to **divide**




Flaherty et al. 2011. Nature Reviews Drug Discovery  
10:811-812

# BRAF-V600E Mutation

- V600 mutation  $\Rightarrow$   
Constitutively active  $\Rightarrow$   
**Uncontrolled cell growth  $\Rightarrow$   
Tumour**
- 50% of melanoma tumours, 10% of colorectal cancers


BRAF codon	599			600			601		
Wild type	A	C	A	G	T	G	A	A	A
V600E	A	C	A	G	A	G	A	A	A
V600K	A	C	A	A	A	G	A	A	A
V600D	A	C	A	G	A	T	A	A	A
V600R	A	C	A	A	G	G	A	A	A
V600G	A	C	A	G	G	G	A	A	A
V600M	A	C	A	A	T	G	A	A	A



Flaherty et al. 2011. Nature Reviews Drug Discovery 10:811-812

# ***BRAF*-V600 Mutation Tests**

- Presence/absence of MT-*BRAF* affects treatment
  - Melanoma patients with mutation can take vemurafenib - a *BRAF* inhibitor



Wagle et al. "Dissecting therapeutic resistance to RAF inhibition in melanoma by tumor genomic profiling"


# ***BRAF-V600* Mutation Tests**

- Cobas® 4800 BRAF V600 Mutation Test (Roche)
- Only looks for V600E
- Requires > 5% mutation level
- Output is qualitative: yes/no





Cobas® 4800  
BRAF  
V600 Mutation Test

# Droplet Digital PCR (ddPCR)



# BRAF-V600 Mutation ddPCR Assay



$$BRAF \text{ mutation frequency} = \frac{\# \text{ MT droplets}}{\# \text{ MT droplets} + \# \text{ WT droplets}}$$



# ***BRAF-V600* Mutation ddPCR Assay**

**Expectation**



**vs**

**Reality**

**Big Mac**





**ddPCR**



# Gating ddPCR Data

- Usually done manually
- QuantaSoft (official analysis software of ddPCR) has auto gating
  - Often wildly inaccurate
- Two tools developed for automatic analysis
  - Both work on single-channel data only
  - Both rely on representative control samples



# Gating ddPCR Data

Beaver, Julia A., et al. "Detection of cancer DNA in plasma of patients with early-stage breast cancer." *Clinical Cancer Research* 20.10 (2014): 2643-2650.

“Droplets were scored as positive or negative based upon their fluorescence intensity which was **determined by gating a threshold** using positive and negative controls as well as no template controls”

Roberts, Chrissy H., et al. "Killer-cell Immunoglobulin-like Receptor gene linkage and copy number variation analysis by droplet digital PCR." *Genome Med* 6 (2014): 20.

“**Crosshair gating was used** to split the data into four quadrants”

Pretto, Dalyir, et al. "Screening Newborn Blood Spots for 22q11. 2 Deletion Syndrome Using Multiplex Droplet Digital PCR." *Clinical chemistry* 61.1 (2015): 182-190.

“The QuantaSoft software (version 1.4.0.99) includes a **freedraw tool** that enables proper classification of the multiple clusters”

Taly, Valerie, et al. "Multiplex picodroplet digital PCR to detect KRAS mutations in circulating DNA from the plasma of colorectal cancer patients." *Clinical chemistry* 59.12 (2013): 1722-1731.

“The sizes and locations of the wild-type gate and the mutant gate(s) were **established by manual selection of the area containing wild-type or mutant clusters**”

Milbury, Coren A., et al. "Determining lower limits of detection of digital PCR assays for cancer-related gene mutations." *Biomolecular Detection and Quantification* 1.1 (2014): 8-22.

“**Objective automated gating** of droplet event clusters is likely **necessary for dPCR practitioners** to take advantage of the full potential sensitivity of the technology for routine applications”


# Goal 1 - Gating ddPCR Automatically

- Given ddPCR output  $\Rightarrow$  calculate mut*BRAF* frequency
- Objective
- Reproducible
- Better gating than QuantaSoft
- No such tools exist

# Goal 2 - Make it easily accessible


- Make R package
  - For people comfortable with R
  - *ddpccr* (on CRAN)
- Make web application with visual UI
  - For people who want a point-n-click interface
  - Uses R package under the hood
  - <http://daattali.com/shiny/ddpccr>

# Dataset: FFPE from 41 CRC Patients




Factors to consider:

- Rain




# Dataset: FFPE from 41 CRC Patients



Factors to consider:


- Failed runs (e.g. C10)

# Dataset: FFPE from 41 CRC Patients




Factors to consider:

- Can't use same thresholds globally (e.g. C08 vs C09)






# Dataset: FFPE from 41 CRC Patients




Factors to consider:

- Outliers




# Dataset: FFPE from 41 CRC Patients



Factors to consider:



- Most droplets are empty

F06;  $\alpha = 0.05$




# Goal

C08




# First try: off-the-shelf clustering algo's


Raw data




K-means




GMM




Hierarchical clustering



Finding two linear equations



My tool




# General Pipeline

1. Identify failed experiments
2. Identify outlier droplets
3. Identify empty droplets
4. Gate droplets (rain vs mutant vs wild type)
5. Classify each sample as mutant or wild type
6. (Revisit gating of wild type samples)





# Step 1: Identify failed experiments

- Use QC metrics to ensure enough data in well





# Step 2: Identify outlier droplets

Take top x% of droplets,  
define threshold as  $Q3 + 5IQR$



# Step 3: Identify empty droplets

Fit two-component Gaussian mixture model to FAM values →  
Lower population is empty droplets, threshold =  $\mu + 5\sigma$






# Step 4: Gate droplets (rain/MT/WT)


First substep: remove the rain

Fit two-component GMM to FAM, threshold =  $\mu - 3\sigma$




# Step 4: Gate droplets (rain/MT/WT)


1. Raw data




2. Remove empty and rain



3. Kernel density estimation (KDE) of HEX values




4. Assuming 2 peaks and 1 trough, use trough as gate




**8.2% mut*BRAF***


# KDE bandwidth selection




Bandwidth  
too small



Bandwidth  
too big





Bandwidth  
just right



# KDE bandwidth selection

Start with low bandwidth, if more than 2 peaks, increase



# Step 5: Classify sample as MT/WT

- Mutation frequency statistically significantly  $> 1\% \Rightarrow$  Mutant
- Use binomial test: What's prob. of observing at least N mutant droplets if the true mutant freq is 1%?
- Example: 500 droplets, 7 mutant.  $H_0$ : freq is  $< 1\%$

Prob observing at least 7 mutants

$$= P(X \geq 7)$$

$$= 1 - P(X < 7)$$

$$= 1 - [P(X=0) + P(X=1) + \dots + P(X=6)]$$

$$= 0.237$$

$> p$ value



Well is classified as WT even though  $\text{mut}BRAF = 1.4\%$

$$P(X=r) \\ = {}_n C_r \cdot p^r \cdot q^{n-r}$$

where  $n$  = total droplets,  
 $r$  = mutant droplets,  
 $p = 0.01$  (1%)

# Results: 41 CRC dataset

- All MT/WT classifications agree with pathologist
- Excellent agreement with manual approach (TOST pvalue:  $1.8 \times 10^{-14}$ )
- 64 seconds on my 3 year old personal laptop



# Results: 41 CRC dataset



```
> plate_meta(myplate)
Source: local data frame [96 x 18]
```

	well (chr)	sample (lgl)	row (chr)	col (int)	used (lgl)	drops (int)	success (lgl)	drops_outlier (int)	drops_empty (int)	drops_non_empty (int)
1	A01	NA	A	1	TRUE	14576	TRUE	0	13884	692
2	A02	NA	A	2	TRUE	15509	TRUE	0	14437	1072
3	A03	NA	A	3	TRUE	16309	TRUE	0	15284	1025
4	A04	NA	A	4	TRUE	14860	TRUE	0	14652	208
5	A05	NA	A	5	TRUE	13879	TRUE	0	13273	606
6	A06	NA	A	6	TRUE	14591	TRUE	0	13893	698
7	A07	NA	A	7	TRUE	13868	TRUE	0	13612	256
8	A08	NA	A	8	TRUE	15280	TRUE	0	14637	643
9	A09	NA	A	9	TRUE	14994	TRUE	0	14118	876
10	A10	NA	A	10	TRUE	14126	TRUE	0	13890	236
..	...	...	...	...	...	...	...	...	...	...

variables not shown: drops\_empty\_fraction (dbl), concentration (int), mutant\_border (int),  
filled\_border (int), significant\_mutant\_cluster (lgl), mutant\_num (int), wildtype\_num  
(int), mutant\_freq (dbl)



# QuantaSoft vs *ddp*cr





# Acknowledgements



Jennifer  
Bryan



Charles  
Haynes



Ryan  
Brinkman



Roza  
Bidshahri



PavLab (Paul Pavlidis)




CIHR Strategic Training Program in  
**BIOINFORMATICS**



THE  
UNIVERSITY OF  
BRITISH  
COLUMBIA

# Summary

1. Identify failed experiments
2. Identify outlier droplets
3. Identify empty droplets
4. Gate droplets (rain vs mutant vs wild type)
5. Classify each sample as mutant or wild type
6. (Revisit gating of wild type samples)




R package: ddpcr

Online: <http://daattali.com/shiny/ddpcr/>


# Step 6: Revisit gating of WT samples

- Wells with few MT droplets don't have enough data to accurately gate
- Look at all mutant samples, we have an idea of where mutant drops are relative to wild type drops



# B-Raf active vs inactive states

- Activation loop (orange) has strong hydrophobic interactions with P-loop (grey)
- These interactions keep the kinase inactive
- Activation loop gets phosphorylated → kinase becomes active
- Valine (V) hydrophobic, glutamic acid (E) is hydrophilic
- V600E → hydrophobic interactions are lost → kinase always active




Wikipedia - BRAF

# FFPE

- Formalin-fixed, paraffin-embedded
- A way to preserve tissue DNA
- Alternative to freezing
  - Less ideal, but doesn't take up as much space and more practical
- Treat sample with formalin solution (which contains formaldehyde) to crosslink the DNA and fix it in place, then put it in paraffin wax
- Formalin causes some degradation, and also causes some C > T mutations


# PCR

Amplify a specific piece of target DNA




# qPCR (real-time)

Monitor fluorescence that gets emitted during amplification to quantify starting amount





Source: Bio-Rad

# Hydrolysis probe





# Copies of target / Droplet $\sim$ Poisson



Example: if **20,000 droplets** and **5000 DNA molecules**, expect 0.25 copies / droplet on average  
This means 78% of droplets will be empty, 19.5% will have one template, virtually none will have 4+

# definetherain

- Upload all positive well
- Use kmeans to define a threshold for positive and threshold for negative (center  $\pm 3$  SD)
- Upload negative wells, and it will use the same thresholds to define positive, negative, and rain
- Concentration is calculated without including rain
- Assumes all wells have very similar distribution
- Still requires manual work of deciding which wells positive and which negative, & upload in two batches
- Kmeans fails if there is lots of rain

# ddpcRquant

- Use combined data of multiple NTCs to model the extreme values of negative droplets by extreme value theory and set a threshold based on that
- Threshold is defined as the 99.5 percentile of the fitted extreme value distribution & used to classify negative threshold in every well
- Droplets are assigned to k groups (blocks) → maximum fluorescence intensity in each group (called the block maxima method) is used to estimate parameters for a generalized extreme value distribution → this distribution used to define threshold
- Assume all wells have same distribution of negatives as NTC
- Website claims R package available Nov 2015, still just an R script

# Poisson to calculate concentration

$$P(x, u) = (e^{-u})(u^x)/x!$$

If we set  $x = 0$ , then

$P(0, u)$  = prob droplet contains no templates

$$= e^{-u}$$

= fraction of negative droplets

$P(x, u)$  = chance of having  $x$  copies in a droplet (where  $x$  = number of copies in a droplet,  $u$  = CPD)

So if we know how many negative/positive droplets we have, we can use poisson equation to figure out the  $u$  (average copies per droplet) in the sample


$$q = e^{-u} \rightarrow -\ln(q) = u \quad (q = \text{fraction of negative droplets})$$

For example, if 75% of droplets are empty, then CPD is  $-\ln(.75) = 0.288$

If the droplet had a total of  $N$  droplets, then  $0.288 * N$  = total copies of target in initial sample




# Step 1: Failed wells conditions

1. # droplets > threshold parameter
2. Empty and non-empty cluster must be well-separated




3. Empty cluster must be not too big nor too small

# Step 4: More examples




# Step 4: Heuristics


If only one peak initially,  
assume all droplets are WT




Every iteration, if <10% of  
droplets are beyond right-  
most peak, discard it



New gate is calculated as  
center+3SD of mutants, if it's  
closer then use it instead





# Results: Horizon samples



Real MT freq	Calculate MT freq
1.4	1.89
0.8	0.691
0	0
50	48.2



# Results: CRC repeat



# Results: brafv600k\_plasmid\_cellline

